

Product datasheet for **MR200106**

Pkig (NM_011106) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pkig (NM_011106) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Pkig
Synonyms: PKIgamma
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200106 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGATGGAAGTCGAGTCCTCTACTCGGACTTCATCTCCTGCGACCGGACAGGCCGTCGGAATGCAGTCC
 CTGACATCCAGGGTGACTCGGAGCCGTGAGTGTGCGGAAGCTCGCTGGAGACATGGGCGAGCTCGCACT
 TGAGGGAGCAGAAGGACAGGCAGAGGGAAGTACCCCGACAAGGAAGCCAGCAGCCAGCCTGAGAGCAGT
 GATCGAACACCTCATCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200106 protein sequence
 Red=Cloning site Green=Tags(s)
 MMEVESSYDFISCDRTGRRNAVPIQGDSEAVSVRKLAGDMGELALEGAEGQAEGSTPDKEASSQPES
 DANTSS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_011106

ORF Size: 228 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_011106.2](#), [NP_035236.1](#)

RefSeq Size: 1066 bp

RefSeq ORF: 231 bp

Locus ID: 18769


UniProt ID: [O70139](#)

Cytogenetics: 2 84.4 cM

MW: 7.9 kDa

Gene Summary: Extremely potent competitive inhibitor of cAMP-dependent protein kinase activity, this protein interacts with the catalytic subunit of the enzyme after the cAMP-induced dissociation of its regulatory chains.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR200106