

## Product datasheet for **MR200108L4V**

### Hsbp1 (NM\_024219) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hsbp1 (NM_024219) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hsbp1
Synonyms:	0610007A03Rik; AA408769; Hsp25
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024219
ORF Size:	231 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200108).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024219.1</a> , <a href="#">NP_077181.1</a>
RefSeq Size:	1178 bp
RefSeq ORF:	231 bp
Locus ID:	68196
UniProt ID:	<a href="#">Q9CQZ1</a>
Cytogenetics:	8 E1
Gene Summary:	Negative regulator of the heat shock response. Negatively affects HSF1 DNA-binding activity. May have a role in the suppression of the activation of the stress response during the aging process (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)