

Product datasheet for **MR200123A1V**

Mouse **BC061195 (BC061195) AAV Particle**

Product data:

Product Type: AAV Particles
Product Name: Mouse BC061195 (BC061195) AAV Particle
Tag: Myc-DDK
Symbol: BC061195
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200123 representing BC061195
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGAAGCAACTAAGAAGACACACTATCTGACACAGGATGTTGTGCCATCAACTTCGGCATCCAAAC
CCAAGACAGTGAAGATTTGCCATGGTCAGCCCTGTTCCAAATGTAAACAAGAAGATAAGAAACCCCAAGGT
GCAGGAGCCTAAGAAGAAGAGCAAGAAAACACTGCCAAGCATGCCTAAAAGAAAGACTTCGAAGAAAAC
ACACCACAAAAACAAAGGACCAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200123 representing BC061195
Red=Cloning site Green=Tags(s)

MVKATKKTHYLTQDVVPSTSASKPKTVKICHGQPCSKCNKKIRNPKVQEPKKKSKKTLPSMPKRKTSKKT
TPPKTKDQ

TRTRPLEQ**KL**ISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: BC061195
ORF Size: 234 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [BC061195.1](#)

RefSeq Size: 746 bp

RefSeq ORF: 236 bp

Locus ID: 69357

Cytogenetics: X C3

MW: 8.8 kDa