

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_183259.3</u> , <u>NP_899082.1</u>
RefSeq Size:	956 bp
RefSeq ORF:	237 bp
Locus ID:	66528
UniProt ID:	<u>Q8BT42</u>
Cytogenetics:	11 E2
MW:	9 kDa