

Product datasheet for **MR200157A1V**

Mouse Ctxn1 (NM_183315) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Ctxn1 (NM_183315) AAV Particle
Tag: Myc-DDK
Symbol: Ctxn1
Synonyms: BC028881; Ctxn
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200157 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGTTCGGCGTGGACACTATCACCGAGCCGCTGCCGCCGTCGACGGGGCCCCCGGTGGGCGCGGGCC
 TGGACGTCGAGCAACGCACGGTGTTCGCCTTCGTGCTCGCTGCTCGTGGTGTGGTTCTGTTGATGGT
 GCGCTGTGTGCGCATACTGCTAGACCCCTACAGCCGCATGCCCGCCTCGTCTGGACCGACCACAAGGAG
 GCGCTCGAACCGGGCAGTTCGACTATGCGTTGGTG

ACGCGTACGCGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200157 protein sequence
 Red=Cloning site Green=Tags(s)

MSSAWTLSPEPLPPSTGPPVAGLDVEQRTVFAFVLCLLVVLVLLMVRVIRILLDPYSRMPASSWTDHKE
 ALERGFQFDYALV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_183315
ORF Size: 249 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_183315.1</u>
RefSeq Size:	1242 bp
RefSeq ORF:	249 bp
Locus ID:	330695
UniProt ID:	<u>Q8K129</u>
Cytogenetics:	8 1.99 cM
MW:	9.1 kDa