

Product datasheet for **MR200173L4V**

Gchfr (NM_177157) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gchfr (NM_177157) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gchfr
Synonyms:	2010323F13Rik; GFRP; P35
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_177157
ORF Size:	252 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200173).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177157.4 , NP_796131.1
RefSeq Size:	629 bp
RefSeq ORF:	255 bp
Locus ID:	320415
UniProt ID:	P99025
Cytogenetics:	2 E5
Gene Summary:	Mediates tetrahydrobiopterin inhibition of GTP cyclohydrolase 1. This inhibition is reversed by L-phenylalanine (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)