

## Product datasheet for **MR200176L4V**

### Cdc42se2 (NM\_178626) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdc42se2 (NM_178626) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cdc42se2
Synonyms:	2810404F18Rik; AA408783; AA536669; SPEC2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178626
ORF Size:	252 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200176).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_178626.3</a> , <a href="#">NP_848741.1</a>
RefSeq Size:	3138 bp
RefSeq ORF:	255 bp
Locus ID:	72729
UniProt ID:	<a href="#">Q8BGH7</a>
Cytogenetics:	11 B1.3



[View online »](#)

**Gene Summary:**

Probably involved in the organization of the actin cytoskeleton by acting downstream of CDC42, inducing actin filament assembly. Alters CDC42-induced cell shape changes. In activated T-cells, may play a role in CDC42-mediated F-actin accumulation at the immunological synapse. May play a role in early contractile events in phagocytosis in macrophages (By similarity).[UniProtKB/Swiss-Prot Function]