

Product datasheet for **MR200190L3V**

Srp9 (NM_012058) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Srp9 (NM_012058) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Srp9
Synonyms:	9kDa
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012058
ORF Size:	261 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200190).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012058.2 , NP_036188.1
RefSeq Size:	1388 bp
RefSeq ORF:	261 bp
Locus ID:	27058
UniProt ID:	P49962
Cytogenetics:	1 H5



[View online »](#)

Gene Summary:

Signal-recognition-particle assembly has a crucial role in targeting secretory proteins to the rough endoplasmic reticulum membrane. SRP9 together with SRP14 and the Alu portion of the SRP RNA, constitutes the elongation arrest domain of SRP. The complex of SRP9 and SRP14 is required for SRP RNA binding.[UniProtKB/Swiss-Prot Function]