

Product datasheet for **MR200196L3V**

Blcap (NM_016916) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Blcap (NM_016916) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Blcap
Synonyms:	AI462828; Bc10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016916
ORF Size:	261 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200196).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016916.3 , NP_058612.3
RefSeq Size:	2099 bp
RefSeq ORF:	264 bp
Locus ID:	53619
UniProt ID:	P62951
Cytogenetics:	2 78.4 cM
Gene Summary:	May regulate cell proliferation and coordinate apoptosis and cell cycle progression via a novel mechanism independent of both p53/TP53 and NF-kappa-B.[UniProtKB/Swiss-Prot Function]



[View online »](#)