

## Product datasheet for **MR200198L3V**

### Pde6g (NM\_012065) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Pde6g (NM_012065) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pde6g
Synonyms:	p3AT; Pdeg
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012065
ORF Size:	264 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200198).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012065.2</a>
RefSeq Size:	919 bp
RefSeq ORF:	264 bp
Locus ID:	18588
UniProt ID:	<a href="#">P09174</a>
Cytogenetics:	11 84.14 cM
Gene Summary:	Participates in processes of transmission and amplification of the visual signal. cGMP-PDEs are the effector molecules in G-protein-mediated phototransduction in vertebrate rods and cones.[UniProtKB/Swiss-Prot Function]



[View online »](#)