

## Product datasheet for **MR200224L3V**

### MIlt11 (NM\_019914) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	MIlt11 (NM_019914) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	MIlt11
Synonyms:	Af1q; AI839562; Zfp692
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019914
ORF Size:	273 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200224).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019914.2</a>
RefSeq Size:	2514 bp
RefSeq ORF:	273 bp
Locus ID:	56772
UniProt ID:	<a href="#">P97783</a>
Cytogenetics:	3 F2.1
Gene Summary:	Cofactor for the transcription factor TCF7. Involved in regulation of lymphoid development by driving multipotent hematopoietic progenitor cells towards a T-cell fate.[UniProtKB/Swiss-Prot Function]



[View online »](#)