

## Product datasheet for **MR200240L3V**

### Dpm3 (NM\_026767) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dpm3 (NM_026767) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dpm3
Synonyms:	1110001H19Rik; AV117483
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026767
ORF Size:	276 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200240).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026767.4</a> , <a href="#">NP_081043.1</a>
RefSeq Size:	506 bp
RefSeq ORF:	279 bp
Locus ID:	68563
UniProt ID:	<a href="#">Q9D1Q4</a>
Cytogenetics:	3 F1
Gene Summary:	Stabilizer subunit of the dolichol-phosphate mannose (DPM) synthase complex; tethers catalytic subunit DPM1 to the ER.[UniProtKB/Swiss-Prot Function]



[View online »](#)