

Product datasheet for **MR200243A1V**

Mouse Snrpe (NM_009227) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Snrpe (NM_009227) AAV Particle
Tag: Myc-DDK
Symbol: Snrpe
Synonyms: AL022645; C76690; SME
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200243 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGTACCGCGGCCAGGGCCAAAAGGTGCAGAAGGTGATGGTGCAGCCCATCAACCTTATCTTCAGAT
ACTTGCAAAATAGATCTCGAATTCAGGTGTGGCTGTATGAACAAGTGAATATGCGGATAGAAGGTTGTAT
TATTGGCTTTGATGAGTACATGAACCTCGTATTAGATGATGCAGAAGAGATTCACTTCTAAAACAAAGTCA
AGAAAACAACCTGGGTCGGATCATGCTAAAAGGAGATAATATTACTCTGCTCCAAGTGTTCCAAC

ACGCGTACGCGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200243 protein sequence
 Red=Cloning site Green=Tags(s)

MAYRGQGQKVKVMQPINLIFRYLQNRSIQVWLYEQVMRIEGCIIGFDEYMNVLDDAEEIHSKTKS
RKQLGRIMLKGDNITLLQSVSN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_009227
ORF Size: 279 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_009227.1</u>
RefSeq Size:	508 bp
RefSeq ORF:	279 bp
Locus ID:	20643
UniProt ID:	<u>P62305</u>
Cytogenetics:	1 E4
MW:	10.8 kDa