

Product datasheet for **MR200246L4V**

lapp (NM_010491) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	lapp (NM_010491) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	lapp
Synonyms:	DAP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010491
ORF Size:	279 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200246).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010491.1
RefSeq Size:	690 bp
RefSeq ORF:	282 bp
Locus ID:	15874
UniProt ID:	P12968
Cytogenetics:	6 73.81 cM
Gene Summary:	Selectively inhibits insulin-stimulated glucose utilization and glycogen deposition in muscle, while not affecting adipocyte glucose metabolism.[UniProtKB/Swiss-Prot Function]



[View online »](#)