

Product datasheet for **MR200262A1V**

Mouse Bglap3 (NM_031368) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Bglap3 (NM_031368) AAV Particle
Tag:	Myc-DDK
Symbol:	Bglap3
Synonyms:	AI461847; Bglap-rs1; mOC-X; OC-X; ORG
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_031368
ORF Size:	288 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_031368.3 , NP_112736.2
RefSeq Size:	611 bp
RefSeq ORF:	288 bp
Locus ID:	12095
UniProt ID:	P54615
Cytogenetics:	3 38.81 cM
MW:	10.4 kDa



[View online »](#)