

Product datasheet for **MR200276L3V**

Apoc2 (NM_009695) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Apoc2 (NM_009695) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Apoc2
Synonyms:	MGC117889
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009695
ORF Size:	294 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200276).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009695.2 , NP_033825.1
RefSeq Size:	727 bp
RefSeq ORF:	294 bp
Locus ID:	11813
Cytogenetics:	7 9.94 cM
Gene Summary:	Component of chylomicrons, very low-density lipoproteins (VLDL), low-density lipoproteins (LDL), and high-density lipoproteins (HDL) in plasma. Plays an important role in lipoprotein metabolism as an activator of lipoprotein lipase.[UniProtKB/Swiss-Prot Function]



[View online »](#)