

Product datasheet for **MR200278A1V**

Mouse **Ccl11 (NM_011330) AAV Particle**

Product data:

Product Type: AAV Particles
Product Name: Mouse Ccl11 (NM_011330) AAV Particle
Tag: Myc-DDK
Symbol: Ccl11
Synonyms: eotaxin; Scya11
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200278 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGAGCTCCACAGCGTTCTATTCTGCTGCTCACGGTCACTTCCTTCACCTCCCAGGTGCTGGCTC
ACCCAGGCTCCATCCCAACTTCTGCTGCTTTATCATGACCAGTAAGAAGATCCCAACACACTACTGAA
GAGCTACAAAAGAATCACCAACAACAGATGCACCCTGAAAGCCATAGTCTTCAAGACCAGGTTGGGCAA
GAGATCTGTGCTGACCCCAAGAAGAAGTGGGTCCAGGATGCCACAAAGCACCTGGACCAAAAACCTCCAA
CTCCAAAACCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200278 protein sequence
Red=Cloning site Green=Tags(s)
MQSSTALLFLLLTVTSFVSQVLAHPGSIPTSCCFIMTSKKIPNTLLKSYKRITNNRCTLKAIIVFKTRLGK
EICADPKKKWVQDATKHLQKLPKP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_011330
ORF Size: 294 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_011330.1</u>
RefSeq Size:	1073 bp
RefSeq ORF:	294 bp
Locus ID:	20292
UniProt ID:	<u>P48298</u>
Cytogenetics:	11 49.84 cM
MW:	10.9 kDa