

Product datasheet for **MR200283A1V**

Mouse Smagp (NM_174992) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Smagp (NM_174992) AAV Particle
Tag: Myc-DDK
Symbol: Smagp
Synonyms: MGC7903
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200283 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAACAACCTACCGGCCACACCTTCTCCTGAAGAATTGATGACCACGCCGGTTTTTCAGGCCCTGAGA
 CCCTGTCCCCACAAGCTGAAGAGGCCAGCACAGCGCTCATTGCAGTTGTTATCACCGTGGTGTTCCTCAC
 CCTGCTCTCGGTGGTACCTTGATCTTCTTTCACTTGTACAAGAACAAGGCAGCTACGTCACCTATGAG
 CCTGCAGAAGGGGAGCCCAGCGCCATCCTCCAGATGGAGACTGACTCAGCCAAGGGCAGAGAGAAGGAAG
 AGTACTTCATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200283 protein sequence
 Red=Cloning site Green=Tags(s)

MNLPATPSPEELMTPVFQAPETLSPQAEESTALIAVVITVVFLTLLSVVTLIFFHLYKNKGSYVTYE
 PAEGEPSAILQMETDSAKGREKEYFI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_174992
ORF Size: 294 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_174992.3</u>
RefSeq Size:	1220 bp
RefSeq ORF:	294 bp
Locus ID:	207818
UniProt ID:	<u>Q99KC7</u>
Cytogenetics:	15 F1
MW:	10.7 kDa