

Product datasheet for **MR200292L4V**

Pyv (NM_145435) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pyv (NM_145435) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pyv
Synonyms:	MGC19143
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145435
ORF Size:	297 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200292).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145435.1 , NP_663410.1
RefSeq Size:	511 bp
RefSeq ORF:	297 bp
Locus ID:	217212
UniProt ID:	Q9EPS2
Cytogenetics:	11 D
Gene Summary:	This gut peptide inhibits exocrine pancreatic secretion, has a vasoconstrictory action and inhibits jejunal and colonic mobility.[UniProtKB/Swiss-Prot Function]



[View online »](#)