

## Product datasheet for **MR200293A1V**

### Mouse Scrg1 (NM\_009136) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Scrg1 (NM\_009136) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Scrg1  
**Synonyms:** AW124307  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200293 representing NM\_009136  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGATGAAATCTGTGGTACTTGTTCATCCTTGGGCTAACTTTGCTGTTAGAAACACAAGCCATGCCTTCAA  
GTCGCCTCTCCTGCTACAGAAAGTTGCTAAAGGATCGCAATTGTCACAACCTCCGGAGGGCAGAGCCGA  
CCTGAAGCTGATAGATGCAAATGTCCAGCATCATTTCTGGGATGGGAAGGGATGCGAGATGATCTGCTAC  
TGCAACTTCAGCGAACTGCTCTGCTGCCAAAAGATGTCTTCTTTGGACCAAGATCTCCTTTGTGATCC  
CCTGCAACAACCAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200293 representing NM\_009136  
Red=Cloning site Green=Tags(s)

MMKSVLVILGLTLLLETQAMPSSRLSCYRKLLKDRNCHNLPEGRADLKLIDANVQHFWGDKGCEMICY  
CNFSELLCCPKDVFFGPKISFVIPCNNH

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_009136  
**ORF Size:** 294 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_009136.3, NP_033162.1</u>
<b>RefSeq Size:</b>	700 bp
<b>RefSeq ORF:</b>	297 bp
<b>Locus ID:</b>	20284
<b>UniProt ID:</b>	<u>O88745</u>
<b>Cytogenetics:</b>	8 29.8 cM
<b>MW:</b>	11.6 kDa