

Product datasheet for MR200299L3

Ndufa2 (NM_010885) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ndufa2 (NM_010885) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Ndufa2
Synonyms:	AV000592; B8; C1-B8; CI-B8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200299).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_010885
ORF Size:	300 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_010885.2
RefSeq Size:	602 bp
RefSeq ORF:	300 bp
Locus ID:	17991
UniProt ID:	Q9CQ75
Cytogenetics:	18 B2
Gene Summary:	This gene encodes a subunit of the NADH-ubiquinone oxidoreductase (complex I) enzyme, which is a large, multimeric protein. It is the first enzyme complex in the mitochondrial electron transport chain and catalyzes the transfer of electrons from NADH to the electron acceptor ubiquinone. The proton gradient created by electron transfer drives the conversion of ADP to ATP. The human ortholog of this gene has been characterized, and its structure and redox potential is reported to be similar to that of thioredoxins. It may be involved in regulating complex I activity or assembly via assistance in redox processes. In humans, mutations in this gene are associated with Leigh syndrome, an early-onset progressive neurodegenerative disorder. A pseudogene of this gene is located on chromosome 5. [provided by RefSeq, May 2013]

Product images:



Circular map for MR200299L3