

## Product datasheet for **MR200305A1V**

### Mouse Smim19 (NM\_001012667) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Smim19 (NM\_001012667) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Smim19  
**Synonyms:** AI316807; BG694986; EST-c39  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200305 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGCTGATGACGGTTCTATTGATTACACAGTCCATGAAGCCTGGAACGAAGCCACCAATGTTTACTTGA  
TAGTTATCCTTGTAGCTTTGGTCTCTTCATGTATGCCAAGAGGAATAAAGGAAAATTATGAGAATATT  
TAGCGTGCCTCCACAGAAGGGATGTTGTCAGAGCCAAGCTTTTATGACACAGTGAGCAGAATTCGCTTA  
AGACAACAAGTGGAGGCGCATCCCGTTTCAAGAAAGTATGAGTATCAACAGCCACAAAGCCAGGCCGACA  
GTGTGCAGCTCTATTGGAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200305 protein sequence  
Red=Cloning site Green=Tags(s)

MADDGSIDYTVHEAWNEATNVYLIVILVSFGLFMYAKRNRKIMRIFSVPPTEGMLSEPSFYDTVSRIRL  
RQQVEAHPVSRKYEYQQPQSQADSVQLSLE

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_001012667  
**ORF Size:** 303 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_001012667.1</u> , <u>NP_001012685.1</u>
<b>RefSeq Size:</b>	1476 bp
<b>RefSeq ORF:</b>	339 bp
<b>Locus ID:</b>	102032
<b>UniProt ID:</b>	<u>Q80ZU4</u>
<b>Cytogenetics:</b>	8 A2
<b>MW:</b>	11.6 kDa