

Product datasheet for **MR200321**

Slc36a3 (BC049682) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Slc36a3 (BC049682) Mouse Tagged ORF Clone

Tag: Myc-DDK

Symbol: Slc36a3

Synonyms: PAT3, tramdorin2, TRAMD2

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Cell Selection: Neomycin

ORF Nucleotide Sequence: >MR200321 representing BC049682
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGGGAATGTGCCGCTGCTTAGAGAAGTTGGCAAGTGTACGCGGAATATGTTCCGGTCGGTCTACAGCCT
CCTCAAAGGCAGCAGCAACAGCAGAAGCAGCAGCAGTACCTCCCAAAGAAGGGCCCTCGCAGGGAAGC
TGATGCGCTAATGGAAGCCCTATGGCAGCAGAGGGGTCTCCACCTGATGCTTGGCATTCTCAGAGGCC
TCCTCTAGGAACCCATCTTTCCAGCTGGAAGGGGTGCCACTACCGTTAATCCAGCCCTGCTTGCTC
ACCAAGGAAGTTTGGTAGTAACCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200321 representing BC049682
Red=Cloning site Green=Tags(s)

MGNVPLLREVGKQQRNMFGRSTASSKGSSNSRSSSTSPKKGPRREADALMEALWQQRGSPDPAWHSSEA
SSRNPIFPAGRGATTRLIQPLLAPPRKFGSNP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: BC049682

ORF Size: 306 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC049682.1</u>
RefSeq Size:	1007 bp
RefSeq ORF:	308 bp
Locus ID:	215332
Cytogenetics:	11 B1.3
MW:	36.9 kDa