

## Product datasheet for **MR200329A1V**

### Mouse Dda1 (BC093523) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Dda1 (BC093523) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Dda1  
**Synonyms:** 1500034J01Rik; 1700095J19Rik; R74921  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200329 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAGACTTTCTGAAAGGCTTGCCCGTCTACAACAAGAGCAACTTCAGCAGGTTCCACGCGGACTCTG  
TGTGCAAGGCCTCGAACCGCCGTCCTCAGTATACCTGCCGCCGAGAGTACCCGTGAGAACAGATCAT  
TGTGACAGAAAAGACGAACATACTTCTGCGGTACCTGCACCAGCAATGGGATAAAAAGAATGCTGCCAAG  
AAGAGAGACCAAGAGCAGGTGGAGGCGGAGGGCGAGAGCTCAGCGCCACCTCGCAAGGTGGCCCGGACCG  
ACAGCCCCGACATGCCGAGGACACC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200329 protein sequence  
Red=Cloning site Green=Tags(s)  
MADFLKGLPVYNKSNFSRFHADSVCASNRRPSVYLPTREYPSEQIIVTEKTNILLRYLHQWDKKNAAK  
KRDQEQVEAEGESSAPPRKVARTDSPDMPEDT

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** BC093523  
**ORF Size:** 306 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">BC093523</a></u> , <u><a href="#">AAH93523</a></u>
<b>RefSeq Size:</b>	1897 bp
<b>RefSeq ORF:</b>	308 bp
<b>Locus ID:</b>	66498
<b>Cytogenetics:</b>	8 B3.3
<b>MW:</b>	11.8 kDa