

Product datasheet for **MR200329L4V**

1500034J01Rik (BC093523) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	1500034J01Rik (BC093523) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	1500034J01Rik
Synonyms:	1500034J01Rik; 1700095J19Rik; R74921
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC093523
ORF Size:	306 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200329).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC093523 , AAH93523
RefSeq Size:	1897 bp
RefSeq ORF:	308 bp
Locus ID:	66498
Cytogenetics:	8 B3.3
Gene Summary:	May be involved in ubiquitination and subsequent proteasomal degradation of target proteins. Component of the DDD-E2 complexes which may provide a platform for interaction with CUL4A and WD repeat proteins (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)