

Product datasheet for **MR200341A1V**

Mouse Pmaip1 (NM_021451) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Pmaip1 (NM_021451) AAV Particle
Tag: Myc-DDK
Symbol: Pmaip1
Synonyms: Noxa
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200341 representing NM_021451
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCCGGGAGAAAGGCGCGTCGGAACGCGCCAGTGAACCCAACGCGGGCAGAGCTACCACCTGAGTTCC
CAGCTCAACTCAGGAAGATCGGAGACAAAGTGTATTGCACGTGGAGTGCACCGGACATAACTGTGGTTCT
GGCGCAGATGCCTGGGAAGTCGCAAAAGAGCAGGATGAGGAGCCCAAGCCCAACCCGGGTGCCAGCAGAC
TTGAAGGACGAGTGTGCTCAACTCCGGAGGATTGGAGACAAAGTGAATTTACGGCAGAACTTCTGAATT
TGATTTCAAGCTCTTCAATTTAGTAACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200341 representing NM_021451
Red=Cloning site Green=Tags(s)

MPGRKARRNAPVNPTRAE L PPEFAAQLRKIGDKVYCTWSAPDITV VLAQMPGKSQKSRMRSPSPTRVPAD
LKDECAQLRRIGDKVNL RQKLLN LISKLFNLVT

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_021451
ORF Size: 309 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_021451.2</u> , <u>NP_067426.1</u>
RefSeq Size:	2654 bp
RefSeq ORF:	312 bp
Locus ID:	58801
UniProt ID:	<u>Q9JM54</u>
Cytogenetics:	18 E1
MW:	11.6 kDa