

## Product datasheet for **MR200343A1V**

### Mouse Tceal9 (NM\_011712) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Tceal9 (NM\_011712) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Tceal9  
**Synonyms:** AL024233; Wbp5  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200343 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAAACCCTGCCAAAAATGGAAGGAAACCTAGAAAAGGAGGATGAGCCAAAGCCTGAGGAAGAGCCAA  
AGCCGGAGGAAAAGCCCGAGGAGGGGCAGGAACCCGAGGAGGAGGAGAAATCCGAAGAACTTTTCGGGA  
GAGGCTGATCCGGTCTCTGCAGGATTTTCAGGAGGATATACACAACAGGCACCTGAGCAGTGATGATCTG  
TTTAGAGATGTGGAGGAGCTAGATGAGATCAGGAAAGTGAGAAACAACTTATAGTGACGCGTTGGAAGG  
CTAACCGAAGCCATCCTTACCCTATTTAATG

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200343 protein sequence  
Red=Cloning site Green=Tags(s)  
MKPCQKMEGNLEKEDEPKPEEPEKPEEKPEEGQEPEEEEEKSEETFRERLIRSLQDFQEDIHNRHLSDDL  
FRDVEELDEIRKVRNKLIVTRWKANRSHPPYLM

**SGP**TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_011712  
**ORF Size:** 315 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_011712.2</a></u> , <u><a href="#">NP_035842.1</a></u>
<b>RefSeq Size:</b>	1018 bp
<b>RefSeq ORF:</b>	315 bp
<b>Locus ID:</b>	22381
<b>UniProt ID:</b>	<u><a href="#">Q9DD24</a></u>
<b>Cytogenetics:</b>	X F1
<b>MW:</b>	12.7 kDa