

Product datasheet for **MR200346A1V**

Mouse Mgp (NM_008597) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Mgp (NM_008597) AAV Particle
Tag: Myc-DDK
Symbol: Mgp
Synonyms: Mg; Mglap
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200346 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGAGCCTGCTCCCTCTGGCCATCCTGGCTGCGCTGGCCGTGGCAACCCTGTGCTACGAATCTCAGC
AAAGCATGGAGTCCTATGAAATCAGTCCCTTCATCAACAGGAGAAATGCCAACACCTTTATGTCCCTCA
GCAGAGGTGGCGAGCTAAAGCCAAAAGAGAGTCCAGGAACGCAACAAGCCTGCCTACGAGATCAACAGA
GAGCCCTGCGATGACTACAAGCTGTGTGAGCGCTACGCCATGGTCTACGGCTACAACGCTGCCTACAACC
GCTACTTCAGGCAGCGCCGAGGAGCCAGATAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200346 protein sequence
Red=Cloning site Green=Tags(s)

MKSLPLAILAALAVATLCYESHESMESYEISPFINRRNANTFMSPQQRWRAKAQKRVQERNKPAYEINR
EACDDYKLCERYAMVYGYNAAYNRYFRQRRGARY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_008597
ORF Size: 315 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_008597.3, NP_032623.1</u>
RefSeq Size:	614 bp
RefSeq ORF:	315 bp
Locus ID:	17313
UniProt ID:	<u>P19788</u>
Cytogenetics:	6 G1
MW:	12.4 kDa