

## Product datasheet for **MR200356A1V**

### Mouse Lbh (NM\_029999) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Lbh (NM\_029999) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Lbh  
**Synonyms:** 1810009F10Rik; 6720416L16Rik  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200356 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGCTGTATATTTCCCATTCCTGCTCTGACTATCTGAGATCGGCTGAGATGACCGAGGTGATGATGA  
 ACGCTCCATCCATGGAAGAGATTGGTCTCAGCCCCGCAAGGATGGCCTTTCCTATCAGATCTTCCGGA  
 CCCATCAGACTTTGACCGTTGCTGCAAGCTGAAGGACCGCCTGCCCTCCATAGTGGTGGAAACCACTGAA  
 GGGGAGGTGAAAAGCGGGGAGCTCCGGTGGCCTCCTGAGGAATTCCTGGTCCAGGAGGATGAGCAAGACA  
 ACTGTGAAGAGACAACAATGAAAAGAAGGACCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAAGTTTAA

**Protein Sequence:** >MR200356 protein sequence  
 Red=Cloning site Green=Tags(s)

MSYFPIHCSYLRSAEMTEVMMNAPSMEIIGLSPRKDGLSYQIFPDPSPDFDRCKLKDRLPSIVVEPTE  
 GEVESGELRWPPEEFLVQEDEQDNCEETTNEKKDQ

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_029999  
**ORF Size:** 318 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_029999.3</u>
<b>RefSeq Size:</b>	3068 bp
<b>RefSeq ORF:</b>	318 bp
<b>Locus ID:</b>	77889
<b>UniProt ID:</b>	<u>Q9CX60</u>
<b>Cytogenetics:</b>	17 E1.3
<b>MW:</b>	12.1 kDa