



Product datasheet for MR200384L3

**Lenti ORF clone of MGC:35756 (Myc-DDK-tagged) - Mouse cDNA clone MGC:35756
 IMAGE:4925089 (Pigyl) (BC025114) Mouse Tagged Lenti ORF Clone**

Product data:

Product Type:	Expression Plasmids
Product Name:	Lenti ORF clone of MGC:35756 (Myc-DDK-tagged) - Mouse cDNA clone MGC:35756 IMAGE:4925089 (Pigyl) (BC025114) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Pigyl
Synonyms:	1810008A14Rik; PIG-Y; Pigyl
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200384).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:




* The last codon before the Stop codon of the ORF.



[View online »](#)

Plasmid Map:



ACCN: BC025114

ORF Size: 323 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: BC025114

RefSeq Size: 664 bp

RefSeq ORF: 323 bp

Locus ID: 66268

Cytogenetics: 9 A3

Gene Summary:

This gene encodes a homolog of a human protein that functions in glycosylphosphatidylinositol biosynthesis. The human protein is expressed from an unusual locus that encodes two distinct proteins in upstream and downstream CDSes; however, in mouse these two proteins are expressed from distinct loci. The product of this locus is highly similar to the protein expressed from the human downstream CDS. A separate mouse locus on chromosome 6 is orthologous to the human locus and encodes a protein similar to the human protein expressed from the upstream CDS. [provided by RefSeq, Jul 2008]