

Product datasheet for **MR200410A1V**

Mouse Rprm (NM_023396) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Rprm (NM_023396) AAV Particle
Tag: Myc-DDK
Symbol: Rprm
Synonyms: 2410012A13Rik; Reprimo
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200410 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAATTCAGTGCTGGGCAACCAGACCGACGTGGCTGGCCTGTTCTGGCCAACAGTAGCGAGGCGCTGG
 AGCGCGCGGTGCGCTGCTGCACCCAGGCGTCGGTGGTGACCGACGATGGCTTCGCCGAGGGAGGCCCCGA
 CGAGCGCAGCCTGTACATCATGCGTGTGGTGCAGATCGCAGTCATGTGTGTGCTCTCGCTCACTGTGGTT
 TTTGGCATCTTCTTCTTGGCTGCAACCTGCTCATCAAGTCCGAAGGCATGATCAACTTCCTAGTGAAGG
 ACCGGAGGCCGTCTAAAGAGGTGGAGGCAGTGGTCGTGGGACCCTAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200410 protein sequence
 Red=Cloning site Green=Tags(s)

MNSVLGNQTDVAGLFLANSSEALERAVRCCTQASVVTDDGFAEGGPDERSLYIMRVVQIIVMCVLSLTVV
 FGIFFLGCNLLIKSEGMINFLVKDRRPSKEVEAVVVGPP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_023396
ORF Size: 330 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_023396.3, NP_075885.1</u>
RefSeq Size:	1460 bp
RefSeq ORF:	330 bp
Locus ID:	67874
UniProt ID:	<u>Q9JJ72</u>
Cytogenetics:	2 C1.1
MW:	11.8 kDa