

Product datasheet for **MR200426L4V**

Gtlf3b (BC034332) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gtlf3b (BC034332) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gtlf3b
Synonyms:	AI256713; F3-2; Gm16515; Gtlf3b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC034332
ORF Size:	330 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200426).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC034332 , AAH34332
RefSeq Size:	3322 bp
RefSeq ORF:	332 bp
Locus ID:	24083
Cytogenetics:	11 B2



[View online »](#)