

## Product datasheet for **MR200441A1V**

### Mouse Fam32a (NM\_026455) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Fam32a (NM\_026455) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Fam32a  
**Synonyms:** 2510049I19Rik; AW822030; OTAG-12; Otag12  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200441 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGAGGCCTACGAGCAAGTCCAGAAGGGCCCCCTGAAGCTGAAAGGCGTGGCAGAGCTCGGCGTGACAA  
 AGCGGAAGAAGAAGAAGAAGACAAAGCAAGGCCAAGATGCTGGAGGCCATGGGGACGAGCAAGAAGAG  
 CGAGGAGGAGAAAAGGCGCTGCCTGGACAAGCGCACGCCGGCGCAGGCAGCCTTCGAAAAGATGCAGGAG  
 AAGCGGCAAATGGAAAGAATCCTGAAGAAGGCATCAAAAACCCACAAACAGAGAGTGGAGGACTTCAACC  
 GACACCTGGATACACTCACTGAGCACTATGACATTCCCAAAGTCAGCTGGACCAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200441 protein sequence  
 Red=Cloning site Green=Tags(s)

MEAYEQVQKGPLKLGVAELGVTKRKKKKDKDKAKMLEAMGTSKKEEERCLDKRTPAQAAFEKMQE  
 KRQMERILKKASKTHKQRVDFNRHLDLTLTEHYDIPKVSWTK

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_026455  
**ORF Size:** 339 bp



[View online »](#)

|                      |   |
|----------------------|---|
| <b>Buffer:</b>       | PBS with 0.001% Pluronic F68  |
| <b>Stability:</b>    | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| <b>RefSeq:</b>       | <u>NM_026455.2</u>  |
| <b>RefSeq Size:</b>  | 1921 bp   |
| <b>RefSeq ORF:</b>   | 339 bp  |
| <b>Locus ID:</b>     | 67922   |
| <b>UniProt ID:</b>   | <u>Q9CR80</u>   |
| <b>Cytogenetics:</b> | 8 B3.3  |
| <b>MW:</b>           | 13.2 kDa  |