

## Product datasheet for MR200442

### Bcl7a (BC017640) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Bcl7a (BC017640) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Bcl7a  
**Synonyms:** 4432415N06Rik; AI448316  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR200442 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGCATCGCC

ATGTCGGGCGAGTTCGGTTCGAGCCGAGACCAGGAGCCGGCCAAAGATGATATCAAGAGGGTCATGGCGG  
 CTATCGAGAAAGTGCACAAATGGGAGAAGAAATGGGTGACCGTTGGCGATACATCCCTACGAATCTACAA  
 GTGGGTCCCTGTGACGGAGCCAAAGTTGATGATAAAAACAAGAACAAGAAGAAAGGCAAGGACGAGAAG  
 TGTGGCTCGGAGGTGACCACTCCAGAGAACAGCTCGTCTCTGGGATGATGGACATGCACGGTGAGCCCC  
 GCCCAGCGCTCCCCCTGCCTGGCCTAAGTTGTGTTTTGTTTCTCTAACAGT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200442 protein sequence  
 Red=Cloning site Green=Tags(s)  
 MSGRSVRAETRSRAKDDIKRVMAAIEKVRKWEKKWVTVGDTSLRIYKWVPVTEPKVDDKNKNKKKGKDEK  
 CGSEVTPENSSSPGMMDMHGEP RPSRSPCLGLSCVFVSLTG

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)

## Cloning Scheme:



ACCN: BC017640

ORF Size: 336 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [BC017640](#), [AAH17640](#)

RefSeq Size: 2135 bp

RefSeq ORF: 338 bp

Locus ID: 77045  
 Cytogenetics: 5 F  
 MW: 12.6 kDa

**Product images:**



Circular map for MR200442