

Product datasheet for **MR200446A1V**

Mouse Scgb2b26 (NM_178308) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Scgb2b26 (NM_178308) AAV Particle
Tag: Myc-DDK
Symbol: Scgb2b26
Synonyms: Abpbg4; Abpbg26; Abpg; C2d; Scgb2b4
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200446 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAGGGGACACTTCTCTGCTGGCCTTGCTGGTGATTGGAGAGCTGGGCTTCCAGACAACGGAAGCAT
 GCCTTTCTTTTCCAGAACTTACGGAGCAATACTACTCTAAGGAGGACCTTCTGCATGGTGACCTTTC
 ACAGTTTTATGCTACTGTAGCTGAAAGGGTGGCCTTTGAAAAAATCCAGGATTGCTTTCGTGAGGAAGGA
 CAAAAACCATAATTCTGAATCCCCAAATTATGTTATCTCTATACTTAAGCCAGAATGCAAGAAATACT
 ATGGCAATGACCTCTTAAAGAAAATCAAGATTTCTTAACCGATCAAATATCCAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200446 protein sequence
 Red=Cloning site Green=Tags(s)

MKGTL L L L L L L L L V I G E L G F Q T T E A C L S F A R T Y G A I L T L R R T F L H G D L S Q F Y A T V A E R V A F E K I Q D C F R E E G
 Q K T I I L N P Q I M L S L Y L S P E C K K Y Y G N D L L K K I Q D F L N Q S N I H

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_178308
ORF Size: 339 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_178308.1</u>
RefSeq Size:	529 bp
RefSeq ORF:	339 bp
Locus ID:	110187
Cytogenetics:	7 19.38 cM
MW:	12.8 kDa