

## Product datasheet for **MR200465A1V**

### Mouse Aurkaip1 (BC011245) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Aurkaip1 (BC011245) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Aurkaip1  
**Synonyms:** Aip, Akip  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200465 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGTTCTGGCGCGCCTGACTTCACGGCTGGCCAGGACTGTTGTTCCCTGGGCAGGTTTCAGTCGTTCT  
 GCCCTGGCTCCGGGTGATCGGCAGCTATGCTTTTCGACCCTTTACAGTTTGCAGCCTGCAAGCCCAAG  
 TCGGGCTGCCTCCCTCCCTGGTAAGAGGACCCAGTTAGAGCTTGAGGAGTTCTGGTCCCCAGGAAGATG  
 GCCATCAGTCTCTAGAGAGCTGGCTGACTGCTCAATACTTATTGCCAGACGGAATGTCGAGGTTCCAG  
 TGACTCTGGCTCCCTCCAGTTCTACGAATGTCCACCTAGGCAGGGGGAGGAGGAAGCCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200465 protein sequence  
 Red=Cloning site Green=Tags(s)

MFLARLTSRLARTVVPWAGFSRSCPGSGVIGSYAFRPLYSLQPASPSRAASLP GKRTQLELEEFVPRKM  
 AISPLESWLTAQYLLPRRNVEVPVTLAPSQFYECPPRQGEETAA

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** BC011245  
**ORF Size:** 342 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>BC011245</u> , <u>AAH11245</u>
<b>RefSeq Size:</b>	1082 bp
<b>RefSeq ORF:</b>	344 bp
<b>Locus ID:</b>	66077
<b>Cytogenetics:</b>	4 87.55 cM
<b>MW:</b>	12.7 kDa