

Product datasheet for **MR200485L3V**

PPP3R1 (BC002267) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PPP3R1 (BC002267) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	PPP3R1
Synonyms:	MCIP1, Cnb1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC002267
ORF Size:	347 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200485).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC002267
RefSeq Size:	802 bp
RefSeq ORF:	347 bp
Locus ID:	19058
Cytogenetics:	11 A2
Gene Summary:	Regulatory subunit of calcineurin, a calcium-dependent, calmodulin stimulated protein phosphatase. Confers calcium sensitivity.[UniProtKB/Swiss-Prot Function]



[View online »](#)