

## Product datasheet for **MR200487L4V**

### Mien1 (NM\_025559) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Mien1 (NM_025559) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mien1
Synonyms:	1810046J19Rik; A1463380; Rdx12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025559
ORF Size:	348 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200487).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_025559.2</a> , <a href="#">NP_079835.1</a>
RefSeq Size:	762 bp
RefSeq ORF:	348 bp
Locus ID:	103742
UniProt ID:	<a href="#">Q9CQ86</a>
Cytogenetics:	11 D
Gene Summary:	Increases cell migration by inducing filopodia formation at the leading edge of migrating cells. Plays a role in regulation of apoptosis, possibly through control of CASP3. May be involved in a redox-related process (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)