

Product datasheet for **MR200489**

Ms4a4b (BC028433) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Ms4a4b (BC028433) Mouse Tagged ORF Clone

Tag: Myc-DDK

Symbol: Ms4a4b

Synonyms: A1463180; Chandra; Ly116

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR200489 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGCAAGGACAGGAACAGACCACCATGGCAGTGGTTCCTGGAGTTGCTGTGCCTCAAAGAATTCTGTTA
 TGACATCACAAATGTGGAATGAGAAGAAAGAGAAATCTTGAAGGGGAACCCAAAGTCCTGGGGTTTT
 ACAAGTGATGATTGCTATCATAACCTCAGCTTAGGAATAATAATTTTGACAACCTTATTTCTGAACTA
 CCCACTTCAGTGATGTTAATGGTCCCAATTTGGGGATCAATAATGTTTCATTGTCTCCGGATCCCTGTCCA
 TTGCAGCAGGAGTGACACCTACAAAATGCCTGGTATGTAATTTGGTATGTTTTGCTATGAAAAGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200489 protein sequence
 Red=Cloning site Green=Tags(s)

MQGQEQTMAVVPGVAVPSKNSVMTSQMWNEKKEKFLKGEPKVLGVLQVMIAIINLSLGIILITLTFSEL
 PTVMLMVPIWGSIMFIVSGSLIAAGVTPTKCLVCNLVCFAMEKE

TRTRPLE**QKLISEEDLA**NDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:


ACCN: BC028433

ORF Size: 348 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC028433](#), [AAH28433](#)

RefSeq Size: 747 bp

RefSeq ORF: 350 bp
 Locus ID: 60361
 Cytogenetics: 19 A
 MW: 12.5 kDa

Product images:



Circular map for MR200489