

Product datasheet for **MR200496L4V**

Hebp2 (BC006898) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hebp2 (BC006898) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hebp2
Synonyms:	SOUL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC006898
ORF Size:	348 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200496).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC006898 , AAH06898
RefSeq Size:	986 bp
RefSeq ORF:	350 bp
Locus ID:	56016
Cytogenetics:	10 A3
Gene Summary:	Can promote mitochondrial permeability transition and facilitate necrotic cell death under different types of stress conditions (By similarity). May have low affinity for heme (PubMed:15518569).[UniProtKB/Swiss-Prot Function]



[View online »](#)