

Product datasheet for **MR200497A1V**

Mouse Zfp593 (NM_024215) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Zfp593 (NM_024215) AAV Particle
Tag: Myc-DDK
Symbol: Zfp593
Synonyms: 3110024A21Rik; AV062409; E130106C14Rik; ZNF593
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200497 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAAGCCAAGAAGCGGCGCCGGACCTGGATGAGATTCATCGCGAGCTGAGGCCGAGGGGCTCCCGC
 GGCCCAAGCCGGAGCGGGACGCAGAGCCGACCCGGACCTGCCAGGGGGCGGCCTGCATCGCTGTCTGGC
 TTGCGCGAGGTACTTCATCGATTCAGCCAACCTGAAGACCCACTCCGATCCAAGACCACAAGAAAAGG
 CTGAAGCAGCTGAGTGTCAACCCTACAGTCAGGAAGAAGCTGAGAGGGCAGCGGCATGGGCTCGTATG
 TTCAACCCAGCGGCTGGCGTGCCACAGAAGTGCCACTGATATCCCTGAGATGGACACGTCCACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200497 protein sequence
 Red=Cloning site Green=Tags(s)

MKAKKRRPDLDEIHRELRPQGLPRPKPERDAEPDPLPGGLHRLCLACARYFIDSANLKHFRSKDHKKR
 LKQLSVEPYSQEAEERAAGMGSYVQPQRLGVPTVSTDIPEMDTST

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_024215
ORF Size: 351 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024215.1</u>
RefSeq Size:	1991 bp
RefSeq ORF:	405 bp
Locus ID:	68040
UniProt ID:	<u>Q9DB42</u>
Cytogenetics:	4 D3
MW:	13.2 kDa