

Product datasheet for **MR200502A1V**

Mouse Ccl6 (NM_009139) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Ccl6 (NM_009139) AAV Particle
Tag: Myc-DDK
Symbol: Ccl6
Synonyms: c10; MRP-1; Scya6
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200502 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGAACTCCAAGACTGCCATTTTCATTCTTTATCCTTGTGGCTGTCTTGGTCCCAGGCTGGCCTCA
TACAAGAAATGGAAAAAGAAGATCGTCGCTATAACCTCCAATAATTCATCAAGGCGTTCAAGACACTTC
TTCAGACTGCTGCTTCTCCTATGCCACACAGATCCCGTGTAAGAAAGATTTATATATTTCCCACTAGT
GGTGGGTGCATCAAGCCGGGCATCATCTTTATCAGCAGGAGGGGAACCCAGGTCTGTGCCACCCAGCG
ATCGGAGAGTTCAAGGTCCTAAGCACCTGAAGCAAGGCCAAGATCTAGGAACAAGGTCATTGCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200502 protein sequence
Red=Cloning site Green=Tags(s)

MRNSKTAISFFILVAVLGSQAGLIQEMEKEKDRRYNPPIIHQGVQDTSSDCCFSYATQIPCKRFIYYFPTS
GGCIKPGIIFISRRGTQVCADPDRRVQRCLSTLKQGPRSRNKVIA

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_009139
ORF Size: 351 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_009139.3</u> , <u>NP_033165.1</u>
RefSeq Size:	1436 bp
RefSeq ORF:	351 bp
Locus ID:	20305
UniProt ID:	<u>P27784</u>
Cytogenetics:	11 50.85 cM
MW:	13 kDa