

Product datasheet for **MR200534**

Vamp1 (NM_009496) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vamp1 (NM_009496) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Vamp1
Synonyms: Iew; Syb-1; Syb1; VAMP-1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200534 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGCTCCAGCTCAACCGCCTGCTGAAGGGACAGAAGGAGCTGCCCTGGTGGGGTCTCCTGGTC
CTCCTCCCAACATGACCAGTAACAGGCGGTTACAGCAGACCCAGGCACAAGTGGAGGAGGTGGTGGACAT
CATGCGTGTGAATGTGGACAAGGTCTTGGAGAGGGACCAGAAGTTGTGAGAGCTGGATGACCGGGCTGAC
GCCTTGCAGGCCGGAGCATCACAATTTGAGAGCAGTGTGCCAAGCTAAAAGGAAGTACTGGTGGAAAA
ACTGCAAGATGATGATCATGCTGGGAGCTATCTGTGCCATCATCGTGGTAGTGATTGTAATCTACTTTT
TACT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200534 protein sequence
Red=Cloning site Green=Tags(s)
MSAPAQPPAEGTEGAAPGGPPGPPNMTSNRRLQQTQAQVEEVVDIMRVNVDKVLERDQKLSLDDRAD
ALQAGASQFESSAAKLKRKYWWKNCKMMIMLGAICAIIVVIVIVYFFT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_009496

ORF Size: 357 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_009496.3](#), [NP_033522.1](#)


RefSeq Size: 2758 bp

RefSeq ORF: 357 bp

Locus ID: 22317

UniProt ID: [Q62442](#)
 Cytogenetics: 6 59.32 cM
 MW: 12.9 kDa
 Gene Summary: Involved in the targeting and/or fusion of transport vesicles to their target membrane.
 [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR200534