

## Product datasheet for MR200535A1V

### Mouse Urah (NM\_029821) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Urah (NM\_029821) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Urah  
**Synonyms:** 1190003J15Rik; 2810420C16Rik; HIUHase  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200535 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCTACCGAGAGCAGTCCCCTGACTACTCACGTTCTAGACACTGCCTCAGGGCTCCCTGCCAAGGCC  
 TCTGCCTCCGGCTGTCCCGCCTGGAGGCCCCCTGCCAGCAGTGGATGGAGCTGAGGACAAGCTACACAAA  
 CCTGGACGGTCTGCTGTCTGGGCTCCTGACACCAAGCCAGATAAAGCCAGGCACCTATAAGCTGTTCTTC  
 GACACAGAGCGCTACTGGAAAGAGCGGGGTCAAGAGAGCTTTTACCCCTATGTAGAGGTGGTTTTCACTA  
 TTACAAAGGAGACCCAGAAGTTCACGTACCTCTGCTGCTGAGCCCATGGTCTACACCACCTACCGGGG  
 GAGT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200535 protein sequence  
 Red=Cloning site Green=Tags(s)

MATESSPLTTHVLDTASGLPAQGLCLRLSRLEAPCQQWMLRRTSYTNLDGRCPGLLTPSQIKPGTYKLF  
 DTERYWKERQESFYPPYEVVFTITKETQKFHVPLLLSPWSYTTYRGS

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** NM\_029821



[View online »](#)

<b>ORF Size:</b>	357 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_029821.1</u>
<b>RefSeq Size:</b>	675 bp
<b>RefSeq ORF:</b>	357 bp
<b>Locus ID:</b>	76974
<b>UniProt ID:</b>	<u>Q9CRB3</u>
<b>Cytogenetics:</b>	7 F4
<b>MW:</b>	13.6 kDa