

Product datasheet for **MR200544A1V**

Mouse Bet1 (NM_009748) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Bet1 (NM_009748) AAV Particle
Tag: Myc-DDK
Symbol: Bet1
Synonyms: AW555236; Bet-1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200544 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGCGTGCAGGCCTGGGTGATGGAGCACCTCCTGGCAGCTATGGAACTATGGCTATGCCAATACTG
GTTATAATGCTTGTGAAGAAGAGAATGACAGACTCACTGAAAGCTGAGAAGCAAAGTAACGGCTATAAA
ATCTCTTTCCATTGAAATAGGCCATGAAGTTAAAAATCAAACAAACTATTAGCTGAAATGGATTCACAA
TTTGATTCTACAACTGGATTTCTAGGTAACCATGGGAAGACTGAAGATTTGTCCAGAGGGAGCCAAA
CAAAATTGCTATGCTATATGATGCTGTTTTTCATTGTTTGTCTTTTTTGTCAATTTATTGGATTATTAAC
TAACTGAGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200544 protein sequence
Red=Cloning site Green=Tags(s)

MRRAGLGDGAPPGSYGNYGyantgynACEEENDRLTESLRskvTAIKSLsIEIGHEVKNQNKLLAEMDSQ
FDSTTGLGKTMGRKILSRGSQTKLLCYMMLFSLFVFFVIYWIILR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_009748



[View online »](#)

| | |
|---------------|---|
| ORF Size: | 357 bp |
| Buffer: | PBS with 0.001% Pluronic F68 |
| Stability: | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| RefSeq: | <u>NM_009748.2</u> |
| RefSeq Size: | 1492 bp |
| RefSeq ORF: | 357 bp |
| Locus ID: | 12068 |
| UniProt ID: | <u>O35623</u> |
| Cytogenetics: | 6 A1 |
| MW: | 13.3 kDa |