

## Product datasheet for **MR200544L4V**

### Bet1 (NM\_009748) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Bet1 (NM_009748) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bet1
Synonyms:	AW555236; Bet-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009748
ORF Size:	357 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200544).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009748.2</a>
RefSeq Size:	1492 bp
RefSeq ORF:	357 bp
Locus ID:	12068
UniProt ID:	<a href="#">O35623</a>
Cytogenetics:	6 A1
Gene Summary:	Required for vesicular transport from the ER to the Golgi complex. Functions as a SNARE involved in the docking process of ER-derived vesicles with the cis-Golgi membrane (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)