

## Product datasheet for **MR200548A1V**

### Mouse Dppa5a (BC092353) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Dppa5a (BC092353) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Dppa5a  
**Synonyms:** AA536857; Dppa5; ecat2; Esg1  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR200548 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGATGGTGACCCTCGTGACCCGTAAGATATCCCCCGTGGGTGAAAGTTCCTGAAGACCTGAAAGATC  
 CAGAAGTATCCAGGTCCAGTCGCTGGTGTGAAATATCTGTTGGCCACAGGGATCTCGAATGTCTCA  
 CATCGAGCAGGTGAGCCAGGCCATGTTTGAGCTGAAGAACCCTGGAATCTCCGAAGAACTATCGAGGTC  
 TTCATTTACGGCTCTCAAAACAACAAGATTCGGGCTAAATGGATGCTTCAGTCCATGGCTGAGAGGTACC  
 ACCTGCCAGCAAAAAGGAGTGCTGAAGCTGGAGGAATCCATGAAGACCCTGGAGCTAGGCCAGTGTAT  
 CGAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR200548 protein sequence  
 Red=Cloning site Green=Tags(s)

MMVTLVTRKDIPPWKVPEDLKDPEVFVQVQLVLYLFGPQGSRMSHIEQVVSQAMFELKNLESPEELIEV  
 FIYGSQNNKIRAKWMLQSMARYHLRQQKGVLLKLEESMKTLELQGCIE

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Mouse  
**Serotype:** AAV-2  
**ACCN:** BC092353



[View online »](#)

<b>ORF Size:</b>	354 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">BC092353</a></u> , <u><a href="#">AAH92353</a></u>
<b>RefSeq Size:</b>	597 bp
<b>RefSeq ORF:</b>	356 bp
<b>Locus ID:</b>	434423
<b>Cytogenetics:</b>	9 E1
<b>MW:</b>	13.8 kDa