

Product datasheet for MR200561

Rps20 (NM_026147) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rps20 (NM_026147) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rps20
Synonyms: 4632426K06Rik; Dsk; Dsk4
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200561 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCATTAAAGATACCGAAAGACGCCCGTGGAGCCCGAAGTGGCGATTCACCGAATTCGAATCACGC
 TCACCAGCCGCAACGTGAAGTCGCTGGAGAAGTTTGTGCGGACTTGATCAGAGGCGCCAAGGAAAAAGAA
 TCTGAAAGTGAAAGGACCGGTGCGCATGCCTACCAAGACTTTGAGAATCACTACCAGAAAAACACCTTGT
 GGTGAAGTTCCAAGACCTGGATCGATTCCAGATGAGGATCCACAAGCGACTCATTGATTTACATAGTC
 CTTCTGAGATTGTTAAGCAGATTACTTCCATCAGTATTGAGCCGGGAGTTGAGGTTGAAGTCACCATTCG
 AGATGCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200561 protein sequence
 Red=Cloning site Green=Tags(s)
 MAFKDTGKTPVEPEVAIHRIRITLTSRNVKSLEKVCADLIRGAKEKNLKVKGPRMPTKTLRITTRKTPC
 GEGSKTWDRFQMRIHKRLIDLHSPSEIVKQITSISIEPGVEVEVTIADA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_026147

ORF Size: 360 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026147.6](#), [NP_080423.1](#)

RefSeq Size: 502 bp

RefSeq ORF: 360 bp

Locus ID: 67427


UniProt ID: [P60867](#)

Cytogenetics: 4 A1

MW: 13.4 kDa

Gene Summary: Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Dec 2015]

Product images:



Circular map for MR200561