

Product datasheet for **MR200566L3V**

Ccl24 (NM_019577) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ccl24 (NM_019577) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ccl24
Synonyms:	CKb-6; MPIF-2; Scya24
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019577
ORF Size:	360 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200566).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019577.3
RefSeq Size:	1289 bp
RefSeq ORF:	360 bp
Locus ID:	56221
UniProt ID:	Q9JKC0
Cytogenetics:	5 G2
Gene Summary:	Chemotactic for resting T-lymphocytes, and eosinophils. Has lower chemotactic activity for neutrophils but none for monocytes and activated lymphocytes. Binds to CCR3 (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)