

Product datasheet for **MR200593**

1700012B09Rik (NM_029306) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700012B09Rik (NM_029306) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: 1700012B09Rik
Synonyms: 1700007P19Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR200593 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGGCGGCAGGAGGCGCTGGCAGGGACAGCCTCAGAGGCGGGCCGGAAGGGGAGCAGCCTAGGCCGG
 CAGGGTGGGGTCCCGCACCCGTGCTGAGCCTGGTGGCCCCAAGAGTCTCGCCAACAATGGAAAACATT
 TCTATACTGTGAACCGCATAAGAGAATCAAGGAGGTTTTGGAAGAAGAACTTCAATCAAGAGAGATGAA
 TGCCATGTGAAAAGCCACCTACAGTGGCCCTGGATGGAATCTGGAGCATTAGAAGGAATCTCCCTGTGG
 GAGGCACCATATCAGGACAGCAGAGCAGAAACAGCTCACTCCACAAGCCAAGTACTACTCGAGGCACGG
 AGGACTGAGAAGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200593 protein sequence
 Red=Cloning site Green=Tags(s)
 MRRQEALAGTASEAGREGQPRPAGLGCRTRAEPGGPQESRQWKTFLYCEPHKRIKEVLEEELSIKRDE
 CHVKSPPTVALDGIWSIRRNLPGVGTISGQQSRNSSLPAKYYSRHGGLRR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_029306

ORF Size: 366 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_029306.3](#), [NP_083582.2](#)

RefSeq Size: 696 bp

RefSeq ORF: 366 bp


Locus ID: 69325

UniProt ID: Q9DAE7

Cytogenetics: 9 A2

MW: 13.6 kDa

Product images:



Circular map for MR200593